

## H1FNT-H1 histone family, member N, testis-specific Gene

Reference: PTXBC118635

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	H1.7; testis-specific H1 histone; HANP1; haploid germ cell-specific nuclear protein 1; histone H1t2; testis secretory sperm-binding protein Li 223n; H1 histone family member N, testis specific		

**For research use only.**

---

<b>GENE DESCRIPTION</b>	H1 histone family, member N, testis-specific
<b>GENE ID</b>	341567
<b>ORIGIN SPECIES</b>	Human
<b>ACCESSIONS</b>	BC118635
<b>VECTOR</b>	pDONR223
<b>NCBI SYMBOL</b>	H1FNT
<b>LENGTH</b>	ProteoGenix bp
<b>SEQUENCE PRIMERS</b>	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

### ORF (Only for reference, product provided without STOP codon )

atggaacaggccttgactggtgagcccaagccggtggccccgagagggcggagtgggccatggctgagggcctggccccagtgccgaatcccaggacactcagccactcagctgccagcggaaaaactgtcggggaccatcgaggggctctcaagctccgtctc  
agagtgcccagttggtgctccaggccatctccactcacaagggctgactctggcagctctcaagaaggagctccgaaacgccggctacgaaftgcgaggaagagcggccgccacgaagcggccagggggcagggccaaggccacgctctcgggtcagcggcagcgacgcc  
gccgctactcagggtctggaaggttcccaagcccaggagaaagccgggacgcgcgaggaagaggaggcagcgcgctccctggaggacccccagccgccccggagctccggaggcggccagcccctcgcaagcggccaggaaggccagagaagtggtggag  
acagaacgcggagggcgaagccaaaggccaatgccagggcgaggagaccaggggcaagggccgagagccaaaggagccgctgtgccaagccaaggaggaagcgggagcagacagcggcagagagggggcaggacagggcgtgaaggaagaccacggcggag  
tcagggaatgacaagagcgaagctccaagcccagggaagagaagcaggccaagaagcccgcacagcggaccatccagtag