

C3orf43-chromosome 3 open reading frame 43 Gene

Reference: PTXBC118636

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	C3orf43; single-pass membrane and coiled-coil domain-containing protein 1; single-pass membrane protein with coiled-coil domains 1		

For research use only.

GENE DESCRIPTION	chromosome 3 open reading frame 43
GENE ID	255798
ORIGIN SPECIES	Human
ACCESSIONS	BC118636
VECTOR	pDONR223
NCBI SYMBOL	C3orf43
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgaacaatgaaccacaacctgatatcctgaaggaggcaatgaaaagagtagaccacaaactccaagcgttagaaacacagtgcaagaactagactcaccaggaataacctgatgcagaaatcgaacatcatagtaaggctttggcaagccaagcagcccaagatgagatgtgg
acagcagttcgggcactccagctcacttcaatggaattgaatattttatcacagctacgtcattgaagtactatctgcttcatactcgtgtgcttgagaagctccagacctgggtgagaggtcttccaaccttagcctctgtactcagaagaaaagttaagaacaagcgcgttagagttgtatggg
agtccatactggaggagtggtggcgtgcaagaaggagacatcacagcactttgtacctcttftatgcacgtggttaacaaggcagaacatactgctaaagtgaaggcagatgtacatcagggatgtcacgttctaattaacaatggttaagaaccaggctctgcaggacagttgtctgag
ggctgtgcaggttaattgagaagggaagcagttaggacccctgaaaagcaaaagtcacccctcgaagagttgataccatctgtcaaaaactaa