

FLJ13224-hypothetical protein FLJ13224 Gene

Reference: PTXBC118637

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	uncharacterized LOC79857		

For research use only.

GENE DESCRIPTION	hypothetical protein FLJ13224
GENE ID	79857
ORIGIN SPECIES	Human
ACCESSIONS	BC118637
VECTOR	pDONR223
NCBI SYMBOL	FLJ13224
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgagaaaccccgagtagtcccccacgggggagagttggcgacttcaactctgttttaaaactggcaagacgccatcatcagcaaacacattgtctgccggccgacgtgccgacggcgcccatccacggctacgcgggctggacattcggacaaggaccag
ggtcgcagcccgaaagcacgtggtctgtacgttagagcacaagcagtgccctgcggagttcactgactcccgccggtacacaaacgcgcgcctccatctccagccaagtgtggctccctggccatttacagcgcggacggcccccataaacgggaagccctgcggtgtcg
cccacactacgtacaggtgaacgcaccgggagcaaaacttcgggtccgaaagcaccggacaccagaggtcttggcgcggtccgaagaggcagggcagtggaagacgggctgttttaa