

LRRC37B2-leucine rich repeat containing 37, member B2 Gene

Reference: PTXBC118647

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms			
For research use only.			

GENE DESCRIPTION	leucine rich repeat containing 37, member B2
GENE ID	147172
ORIGIN SPECIES	Human
ACCESSIONS	BC118647
VECTOR	pDONR223
NCBI SYMBOL	LRRC37B2
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgaagatgttacaagcccggaagcagcacatgagcactgagctgactattgagtcggaggccctcagacagcagtgccatcaactgtcaggcttgggagtgagcagctagacacaaatgacgagagtgatgttagcagcgactaagttacatcttgccttatctctcactgagaaa
tctaggatgacagcttgaattcagtaaccgagcagctacggctccctatccccaacgaggatgtgagaaagtcatgtctcatgttatctggaccttgaataatgttcagaacacatgtgcaaggagagctgtgccaagctcatgtcgcgaacaggcctcctgatgaagctctcagc
gagcagcagggaagcaaggcattgaatgtagaattggatagcggaccaacaaaaacaaattatattaatgagaacatggaacagaatgaacagaagagcagaagtcgaatgagctcatgaaagaagttccaggatagactataagaacaaactatcttcgaatatctgtgactgtca
tactaataatttgattataattttgttttatagaggtaaagacaataaattcagggtttcaaaatacaatcctgtgtgtgtggtgattcagaatccacaactgaaaaccaacgtcacttcccacttgacattcttctgtcatttaaggctgaggtgtgtgtgttttactgcaatgtatattcca
ggattgttaa