

C20orf69-chromosome 20 open reading frame 69 Gene

Reference: PTXBC118988

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	C20orf69; NCRNA00266; NCRNA00266-1; bA476I15.3; long intergenic non-protein coding RNA 266-1		

For research use only.

GENE DESCRIPTION	chromosome 20 open reading frame 69
GENE ID	140849
ORIGIN SPECIES	Human
ACCESSIONS	BC118988
VECTOR	pDONR223
NCBI SYMBOL	C20orf69
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgattcaacaggaggagataaggagctcgaggagagagaaaaacaactggaggagaaatcatagattttataaaatgaaagctgcctctgaagcactgcagactcagctgagcactgatacaaaagagaacacatcctgatccatgaattcctcttattaagaaaaataagca
tccaggattcaatgaagaactatcacctgttaa