

## GFOD1-glucose-fructose oxidoreductase domain containing 1 Gene

Reference: PTXBC119005

|                |                                                                                                                                    |                     |                |
|----------------|------------------------------------------------------------------------------------------------------------------------------------|---------------------|----------------|
| Product type:  | DNA                                                                                                                                | Storage Condition:  | -20°C          |
| Size/Quantity: | 2 ug                                                                                                                               | Delivery Condition: | +4°C, Blue Ice |
| Synonyms       | ADG-90; C6orf114; glucose-fructose oxidoreductase domain-containing protein 1; glucose-fructose oxidoreductase domain containing 1 |                     |                |

**For research use only.**

---

|                         |                                                                                                  |
|-------------------------|--------------------------------------------------------------------------------------------------|
| <b>GENE DESCRIPTION</b> | glucose-fructose oxidoreductase domain containing 1                                              |
| <b>GENE ID</b>          | 54438                                                                                            |
| <b>ORIGIN SPECIES</b>   | Human                                                                                            |
| <b>ACCESSIONS</b>       | BC119005                                                                                         |
| <b>VECTOR</b>           | pDONR223                                                                                         |
| <b>NCBI SYMBOL</b>      | GFOD1                                                                                            |
| <b>LENGTH</b>           | ProteoGenix bp                                                                                   |
| <b>SEQUENCE PRIMERS</b> | Forward primer M13R (5'→3'): GTAAACGACGGCCAGT<br>Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG |

### ORF (Only for reference, product provided without STOP codon)

atgttccgggggtggcggttgcggcaccagctcagcgccgtgtcatccgctgtgaaagacgagggcttcgctggaagcgctgtggggccgcacgcagggaagcggaggagctggccaaggagatgagtgtcccttctacactagccgcatgaggtgctg  
ctgcatcaggagctggactgtgtgctcattaacctgcgcggccctcaccagacagatcgctgtcaaaaccctaggcatcggaagaacgtcatctgcgacgcacggccacgcctgtggacgtttccgcatgacctcggcccccactactacccaagctcatgagcatatggg  
caacgtgtcgtcgttctcgtcggcgttctgtgctgaagcagctgatcagaggagggctacgtggcgagccgctgtgtgtgaggtgacgtgacggcgagcctgtggcgaagaagtacaactggagctgcgacgactgatggcgggcgccgctgacatccgtgggca  
cctacatcatcgactgtcctcacttctcaccggccaaaggccgtcaaggtccacgggctgtcaagaccttctggaagcagactgaccacatcaaggcgtatccgacagatcaccagcgtgactctgcaccttcagatggtgctggagggcggggtgtgtgctgacgtcacctca  
actcaacgtgcccggcgagttcaagcaggatgtcactgtgtgtggctcagccggcgctgtgctgcccgtgggcacccgactgtacgggcagcgcaacagcggccggagcaggagctgtgtgtagggagccacggcggtgagcaactcctgttccggagaagccttcag  
cgacatccctcgcctacgtcgcggcaccatcaagatgatgcagcgggtgcgccagcccttcaggaccaggacgaccggcgacgtggatggcgccctcaccatggccacacttcagcagctgctgtatgccttgtgctgtgtgacacatcaagaggtccagcca  
gacggcgagtgccagaacattgccatcatgaccgaggagccggagctgagcccgccctacgtgatcagcgagccatgcggcgacgaggtatcccttactgttag