

C15orf53-chromosome 15 open reading frame 53 Gene

Reference: PTXBC119016

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	uncharacterized protein C15orf53; chromosome 15 open reading frame 53		
For research use only.			

GENE DESCRIPTION	chromosome 15 open reading frame 53
GENE ID	400359
ORIGIN SPECIES	Human
ACCESSIONS	BC119016
VECTOR	pDONR223
NCBI SYMBOL	C15orf53
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atggagctacaagggcccaagaggacctgggcatctccctctctagtcctccggaggaaacatgaacaggccagggaagcaaggcctaaggcaggagcagcatctgtctcaggcgtctgtttggatggctggaggaaagctgaggtcagagcaagtgaacatcaaccagggc
caccaacaagagctcagggatggaaatctgggagaagatgctccctctctctcctcaagagtcctttgggctggcaggtgatacaagctccagctcatgtctcagacaatgctgggtccaggaaatgcatggatctctgattacaagtttgactccaagagaagcccagaagtgg
cttcagccccctctacctaactgtgccccgccgttcaccactctctctctcgcagccagtgacagatgtgtctgtggggctgctacttgggaaagtcaccaggaggagagcatgccaatctctctctcagatcctctgggggttactttccccacccaactaggccttga