

LOC728276-hypothetical LOC728276 Gene

Reference: PTXBC119018

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	C-type lectin domain family 19 member A		

For research use only.

GENE DESCRIPTION	hypothetical LOC728276
GENE ID	728276
ORIGIN SPECIES	Human
ACCESSIONS	BC119018
VECTOR	pDONR223
NCBI SYMBOL	LOC728276
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgcaagggtggacactgtgggctgcagccttctgacccctccagctgcacaggcccttcacaaacagacatcagtcagtcagccctgccagagctgccccctgcccctgtgccccctgttctggatggagttcaaggccactgtatcgtattctccctctcaataagacctggg
ccgaggccgacctctactgttctgagttctctgtggcagggaagtcgcccaagctggcctcatccacagctggaggagaatgtctttgtatgatccctctgtaacagctgtgtccggcatcccagctgacgtctggacaggccctcatgatcacagacaggtgagaaagcagtgccca
ttgggcccccttggagctccagcaggattctattttgatttaa