

ZNF321-zinc finger protein 321 Gene

Reference: PTXBC119627

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms			
For research use only.			

GENE DESCRIPTION	zinc finger protein 321
GENE ID	399669
ORIGIN SPECIES	Human
ACCESSIONS	BC119627
VECTOR	pDONR223
NCBI SYMBOL	ZNF321
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgatgaaggagttctccacagcgcaaggcaatacagaagtgaacacacagggacattgcaagacatgaaagtcacacattagagattttgctccaggaaattgagaaagatattcataactttgagtttcagtggaagaagaggaaaggaatggtcacgaagcaccatgac
agaaatcaagagttgactggtagtagacagccgacatgatcaaggcatgctggaacaagcctattaaagatcagcttgatccagctttcattgcctgcaactccacatatttcagcctgaatggaaaattggaatcaagttgagaagtcacatcaatgcccttaatttgacat
cccaagaatttctgtagtcccaaaaccgtatttctaataactatgggaataatccctcattcttcattaccacatacaaaaattgggaagtacacatgagagaaaaatctttccaatgtaa