

LRRC3-leucine rich repeat containing 3 Gene

Reference: PTXBC119648

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	C21orf102; leucine-rich repeat-containing protein 3; leucine rich repeat containing 3		
For research use only.			

GENE DESCRIPTION	leucine rich repeat containing 3
GENE ID	81543
ORIGIN SPECIES	Human
ACCESSIONS	BC119648
VECTOR	pDONR223
NCBI SYMBOL	LRRC3
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgggcaaccgtgcgccacctcgcctcctcgtcgtgtgtccacccggagctctgtctctctctctctctcctgcacctggcgccgctgccacagccctgccgggtgccctgaccacgcaggggctgtgctgtctctcagcttgcggggccttcaggaggtccccgagg
acatccccggccaacacgtgctcctgaagctcgtatgccaagaatctccacctccggaggggcttcagcagctgcacggctcagggagctggatctgtctcacaacgccatcagggccatcggtccgccacttcggggcctggccggggcctgcggctcgtggacc
tgttttacaaccgcatccaggagatccccaaaggacgcccgggcaaaatcagcgcgaagatagcctgtccacaacccccctgcactgcagtgccctgcaggagggcctgtgggagctgaagctggaccccactctgtggacgagatcgccacacactcagtcaggagg
agtttgggggaagcctctgttcaggctctggtatcggtggcagcctctgcagcgtccccacaggaccacagacgtggccatgtgtgtcaccatgttcggctgttcgccatgtgtgatcgctacgtcgtgtactatgtgcgccacaaccaggagatgccggaggcacctggagt
acctgaagctctgccagcggcccgctccaaggaccccatcgccggggccctag