

C19orf39-chromosome 19 open reading frame 39 Gene

Reference: PTXBC119677

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	C19orf39; SWS1AP1; ZSWIM7AP1; ATPase SWSAP1; SWS1-associated protein 1; ZSWIM7-associated protein 1; zinc finger, SWIM-type containing 7 associated protein 1; SWIM-type zinc finger 7 associated protein 1		

For research use only.

GENE DESCRIPTION	chromosome 19 open reading frame 39
GENE ID	126074
ORIGIN SPECIES	Human
ACCESSIONS	BC119677
VECTOR	pDONR223
NCBI SYMBOL	C19orf39
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgctgcccggaccgctttgctgctggtacaccagatctggaaaaacagcgctgtatttgcggccctagaggcggcggggaggcccaaggccagtcctctctgacacgaaggcctctcaagcatgccccggggaccgggaacgactctagaccaatgc
gactccagaagatccgttcagfaccaccctcaaccgagagctttccggctcctgtgctctgcccagaggccccggggccagccccctcctctctgctgctcagcggcctagaagagtagctagcgggaagaccagagccccaggaagccgctactcattgccttactctaga
cacagctgcccacttcagccaccggcttggcctggccgggaltgtgggctcatggtggccctccagaccaggagggaagcaggtagtggggacgtcctgcacctggcactgctccagcggtattttctgcccagtgctggtcagccagatgcaccagggtccaggagacacgg
cctccgagcctgctgagccaggcgggctgggccccagaacagatggtgggtgactttccgatcagatggagaatgatgcctcgtggccaccaggctgtgacccagctcaggcaagggttcaagctctggaggccagccctga