

BPESC1-blepharophimosis, epicanthus inversus and ptosis, candidate 1 Gene

Reference: PTXBC119686

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	NCRNA00187; blepharophimosis, epicanthus inversus and ptosis, candidate 1 (non-protein coding)		

For research use only.

GENE DESCRIPTION	blepharophimosis, epicanthus inversus and ptosis, candidate 1
GENE ID	60467
ORIGIN SPECIES	Human
ACCESSIONS	BC119686
VECTOR	pDONR223
NCBI SYMBOL	BPESC1
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgggaggaagagatcaaatcacaatctgggagtgaggagagacacatcctggggtactcaagtccccatcctctgcatcctcttaccctttctaccagtgattggaggagtgctggctggataccatcagacaacgagaacacgattgcaggactgtgtggagg
tatccaggcctgaggccctgtccagagctccccctcactctgtgctggaacaaaatctctctgtgtggccttggctccctagcaggatcccaaacatggatctgccatgctgctcagctctgggtgattctttagctccagctctctga