

MGC2752-hypothetical LOC65996 Gene

Reference: PTXBC119702

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	FLJ18389; MGC111383; MGC149628; hypothetical LOC65996		

For research use only.

GENE DESCRIPTION	hypothetical LOC65996
GENE ID	65996
ORIGIN SPECIES	Human
ACCESSIONS	BC119702
VECTOR	pDONR223
NCBI SYMBOL	MGC2752
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgagtcggccctgggtcaaaagccgagagactcctgacccattgctcgtggcaacgcagccggcgactcaagctgaaggccctgctgggtacccctggagaacccgtgttctcaagggtccttcaagcccaacctgcccttggtgtggtcgcacagaagcct
gggtacagcctgggtgacagagcggagactccgtctcaaaaaaaaaaaaaagcctgggttaacatgagcctctccaggagtggttctgcacttctgctctacagtgaggcgctgcacccagtacgacctccctacaagggtgccccatcctgacagcactcccgccaccc
tgggtgttgagacagagctgctctgctgcgcagcctggagtcgggtggcgcatctcggctcactgcaacctccgctccgggttcacgccattctcctcagcctccaggtagctgggactacaggcggccgccaccagcccggttaatttttttagtagagacgg
ggtttcacggtgttagccagaatggtcgtctcctgacccatga