

LCE2A-late cornified envelope 2A Gene

Reference: PTXBC119707

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	LEP9; late cornified envelope protein 2A; late envelope protein 9; late cornified envelope 2A		

For research use only.

GENE DESCRIPTION	late cornified envelope 2A
GENE ID	353139
ORIGIN SPECIES	Human
ACCESSIONS	BC119707
VECTOR	pDONR223
NCBI SYMBOL	LCE2A
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgtctgccagcaaaaccagcagcagtgccagcccccctcccaagtgcgcccaaaatgccacccaagtgtctccaaagtgcgacctcagtgccagcccccagccacctccagtcctctctctgtgtgtgtccagctctgggggctgctgcgggtccagctctgggggctgctgc
agctctgggggtggcgctgctgcctgagccaccacagggccccgtctctccacggcaccggcaccagagccccgatgtgtgtgagtgtaacctctgggggctgctgctgctgccacagctctggggactgctgctga