



Antibodies and Proteins

ProteoGenix SAS

Espace Européen de l'Entreprise

15 rue de La Haye

67300 Schiltigheim, France

Tel. +33 (0)3 90 20 54 70 - Fax +33 (0)3 88 56 16 88

contact@proteogenix.fr

www.ProteoGenix.science

LY6E-lymphocyte antigen 6 complex, locus E Gene

Reference: PTXBC119708

Product type:

DNA

Storage Condition:

-20°C

Size/Quantity:

2 ug

Delivery Condition:

+4°C, Blue Ice

Synonyms

RIG-E; RIGE; SCA-2; SCA2; TSA-1; lymphocyte antigen 6E; ly-6E; retinoic acid induced gene E; retinoic acid-induced gene E protein; stem cell antigen 2; thymic shared antigen 1; lymphocyte antigen 6 complex, locus E

For research use only.

GENE DESCRIPTION	lymphocyte antigen 6 complex, locus E
GENE ID	4061
ORIGIN SPECIES	Human
ACCESSIONS	BC119708
VECTOR	pDONR223
NCBI SYMBOL	LY6E
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgaagaatcttcttccagtcgtcgtcgtcgtccctctctgggttgagcagagccagctcgtgatgctctctctcttgaaccagaagacaatctgtactcgtgaagccgaccatctgtccgaccagacaactactcgtgctgctgtctagtccggcattgggaatctcgtga
catttggccacagctctgagacaagatctctccggcctccccatccagaaggcgtcaattgtgggttggttccatgggcatcagctcgtccagagcttctgtgcaatttcagtcggccgatggcggcgtcggggcaagctcacctcgtcgtgggtcgggcgtcgtcgtgagctt
ctgcgcccccctctcgtcttctgcccctta