

C17orf91-chromosome 17 open reading frame 91 Gene

Reference: PTXBC119720

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	C17orf91; MIR22 host gene (non-protein coding)		

For research use only.

GENE DESCRIPTION	chromosome 17 open reading frame 91
GENE ID	84981
ORIGIN SPECIES	Human
ACCESSIONS	BC119720
VECTOR	pDONR223
NCBI SYMBOL	C17orf91
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgggctgggaaggtccgaacagcagggtggatgatacgttttggcgaagttggagagcctctgccagattggccagcaaggagcgggttttagattagagacactggctggattgaggagtagaaggctcaacaaccaaggctttctgcctccgagatgtggcaccatag