

CORT-cortistatin Gene

Reference: PTXBC119724

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	CST-14; CST-17; CST-29; cortistatin; cortistatin-14; cortistatin-17; cortistatin-29; prepro-cortistatin; preprocortistatin		

For research use only.

GENE DESCRIPTION	cortistatin
GENE ID	1325
ORIGIN SPECIES	Human
ACCESSIONS	BC119724
VECTOR	pDONR223
NCBI SYMBOL	CORT
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgtatagacacaaaacagctggagattggccttaataaccaccaagctccaaagaagagacccaagtcaccaaaacattgattcagggtgccaggaaggaagagcagcagcagggtgggagagaagctccagtcagcccaagaatgccattgtccccggcctcctgctgc
tgcctgtctccggggccacggccaccgctgccttgccttggagggtggccccacggccgagacagcgagcatatgcagggaagcggcagggaataaggaagcagcctcctgactttctcgttgggtgttggagtgacccagccagtcggggccctcataggagag
gaagccggggagggtggccaggcggcaggaaggcgcacccccagcaatctgcgcggggacagaatgccctgcagggaacttcttctggaagaccttctcctcctgcaataaa