



Antibodies and Proteins

ProteoGenix SAS

Espace Européen de l'Entreprise

15 rue de La Haye

67300 Schiltigheim, France

Tel. +33 (0)3 90 20 54 70 - Fax +33 (0)3 88 56 16 88

contact@proteogenix.fr

www.ProteoGenix.science

FKBP1A-FK506 binding protein 1A, 12kDa Gene

Reference: PTXBC119732

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	PPIase FKBP1A; peptidyl-prolyl cis-trans isomerase FKBP1A; FKBP-12; FKBP-1A; FKBP1; PKC12; PKCI2; PPIASE; 12 kDa FK506-binding protein; 12 kDa FKBP; FK506 binding protein 1A, 12kDa; FK506 binding protein12; FK506-binding protein 1; FK506-binding protein 1		

For research use only.

GENE DESCRIPTION	FK506 binding protein 1A, 12kDa
GENE ID	2280
ORIGIN SPECIES	Human
ACCESSIONS	BC119732
VECTOR	pDONR223
NCBI SYMBOL	FKBP1A
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgggcaggcaacgcgcgtgagggcactaggcagagccgtggaaacccgccagagtcgctgtgtgtgccaccgcgcgcctgcgcgcccccgcgcctcagtcagtcgcgcgcccccatgggagtcaggttgaaaccattccccagagagacggcgccaccttccccaaagcgcggccacgcacctgcgtgtgtgactacacccgggatgtgtgaagatggcaagaagaaatgtatctcccccagcagaacacaccccttaagtattatctagggcaagcagagagtgatccgaggtcggtggaaaggggtttccacgatgatgtgggttcagagagccaaactgactatctccagattatgccctatggctcactgcgcacccagcgcacatccccacacatgaactctctcttcgcatgtggagctctctaaacgtgaatga