

CSAG3-CSAG family, member 3 Gene

Reference: PTXBC119735

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	CSAG3A; CT24.2; chondrosarcoma-associated gene 2/3 protein; CSAG family, member 3A; cancer/testis antigen 24.2; taxol resistance associated gene 3; taxol-resistant-associated gene 3 protein; CSAG family member 3		

For research use only.

GENE DESCRIPTION	CSAG family, member 3
GENE ID	389903
ORIGIN SPECIES	Human
ACCESSIONS	BC119735
VECTOR	pDONR223
NCBI SYMBOL	CSAG3
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgtggatgggctcatcaattagttgaaggtgttaagagaaaagaccaaagggttcctggaaaagggaattctaccacaagactaacataaaatgcgctgtgagttctagcctgctggcctgccttactgtcctggggaggcctggagagaccaggtggactggagtagactgttgaga
gacgctggctgtgtaagatgtccaggaaccacgagcctccagccattgtccaacaaccacccaccaacaccaaagaggttccaagacaacccgaagggaagggaacccatcaaggaggttcaggaaacaaaggctctcctaa