

C21orf88-chromosome 21 open reading frame 88 Gene

Reference: PTXBC119737

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	chromosome 21 open reading frame 88		
For research use only.			

GENE DESCRIPTION	chromosome 21 open reading frame 88
GENE ID	114041
ORIGIN SPECIES	Human
ACCESSIONS	BC119737
VECTOR	pDONR223
NCBI SYMBOL	C21orf88
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgagacgactccggcacagagaagtcagagggcctgtcctgggtcacacagccaggcggccccagaggggactcagggtgcaccactgcacccagcagagacctccgctggcacaggggaatgctggagcaatactaaatagggagggtcagaaatctcacgggc
tctgctggcttctgtgcttcgtctcctcaaggccaaaatcaagatgcatcagtgctgagctctgggtgaggatccacgtccaggctcactgggctgctggcagaatcagctgtgtgggtgttaggaatgaagtctgtttccttgcctgtgtggccaggggctgccctcagcttctggag
gccgctgccctccctggttcattggcccatctcctcagagccagcaccgcgggaagtccctcttgcactgtga