

RNF5-ring finger protein 5 Gene

Reference: PTXBC119741

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	E3 ubiquitin-protein ligase RNF5; RING5; RMA1; protein G16; ram1 homolog; ring finger protein 5, E3 ubiquitin protein ligase; ring finger protein 5		

For research use only.

GENE DESCRIPTION	ring finger protein 5
GENE ID	6048
ORIGIN SPECIES	Human
ACCESSIONS	BC119741
VECTOR	pDONR223
NCBI SYMBOL	RNF5
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atggcagcagcggaggaggagggggggccccaaggggccaaatcgcgagcggggcggggcggcgcgaccttcgaatgtaatatatgtttggagactgctcgggaagctgtgtcagtggtgtgtgtgccacctgtactgttgccatgtcttcacagtggtggtggagacacggcca
gaacggcaagagtggtccagtatgtaaaagctgggatcagcagagagaaggtgtcccgctttatggcgaggaggagccagaagcccgagatccagattaaaaactccacccggcccccaggccagagaccagctccggagagcagaggggattccagccatttggataccgg
gggctccacttcacatttgggtgtgtgtgtcttcccttggcttttaccaccgctctcaatgcccatgagccttccgcccgggtacaggtgtgtgacgtgggacaggtgacccggcctccagctggcaggatccctctctctgtttctcgccatcttcttttttggctgtcagctattga