



Antibodies and Proteins

ProteoGenix SAS

Espace Européen de l'Entreprise

15 rue de La Haye

67300 Schiltigheim, France

Tel. +33 (0)3 90 20 54 70 - Fax +33 (0)3 88 56 16 88

contact@proteogenix.fr

www.ProteoGenix.science

RNF5-ring finger protein 5 Gene

Reference: PTXBC119742

Product type:

DNA

Storage Condition:

-20°C

Size/Quantity:

2 ug

Delivery Condition:

+4°C, Blue Ice

Synonyms

E3 ubiquitin-protein ligase RNF5; RING5; RMA1; protein G16; ram1 homolog; ring finger protein 5, E3 ubiquitin protein ligase; ring finger protein 5

For research use only.

GENE DESCRIPTION	ring finger protein 5
GENE ID	6048
ORIGIN SPECIES	Human
ACCESSIONS	BC119742
VECTOR	pDONR223
NCBI SYMBOL	RNF5
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgcagccagcgggagggagggacggggggcccccgaaggcccaatcgcgagcggggcgggggcgggcggaactcgaatgtaatatgtttggagactcctcgggaagctggtggcagtggtgtgtggccaccctgtactgttggccattcttcatcagtgctggagacacggccac
gaacggccaagatggtccagatgtaaaagctcagacagagagaaagttgtcccgctttatggcgagagcagaagacccagatccagatgaaataaaactcaccgcgcccccggggcgagagactcagagagcagagagggatcttcagcgaatttgcgtgatccgg
ggcttcctactcttcattgctgtgtgtgtcttccctcttgcttttccactcagccttcctccccaatgaccttccgcggcggtacagctgttgatctcgcagatgctcaccgcggccctcagctggccagagatccctctctctcttctctccctcctcttcctgatttgc