

AKIRIN1-akirin 1 Gene

Reference: PTXBC119745

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	C1orf108; STRF2; akirin-1; akirin 1		
For research use only.			

GENE DESCRIPTION	akirin 1
GENE ID	79647
ORIGIN SPECIES	Human
ACCESSIONS	BC119745
VECTOR	pDONR223
NCBI SYMBOL	AKIRIN1
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atggcgtgcggggcgagcgtgaagcggcccatggagttcgaaggcggcgtgctgagcccggtcccgaaagcggggcgtgcgcccctctgccgggccccactccggcctcagccccggagccgagccgcccgtttcagacgcagacccaccgcagagctc
gcagcagccccggccggcagcggcggcgttcacactccggagcaaattttcagaacataaaacaagaatagtcgtatcagaggtggagacattgaagttgttttaacagagtgaaagctgtgtctcggaaagtaaacctcactcctcagcactcacgacactagctc
tccagggttctcatggaagaaggaccagccacatttaccctccgacaagttggcataatgtgagcgcctcttaaaagactatgaagataaaattcgggaggagtgagcaaatcctcaataccaactagcagaacaatgaatctttgtgaaattcacacatgatcagattatgc
gacggtatgggacaaggccaacaagctatgtgtcatga