

CXXC4-CXXC finger 4 Gene

Reference: PTXBC119751

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	IDAX; CXXC-type zinc finger protein 4; CXXC finger 4; Dvl-binding protein IDAX (inhibition of the Dvl and Axin complex); inhibition of the Dvl and axin complex protein; CXXC finger protein 4		

For research use only.

GENE DESCRIPTION	CXXC finger 4
GENE ID	80319
ORIGIN SPECIES	Human
ACCESSIONS	BC119751
VECTOR	pDONR223
NCBI SYMBOL	CXXC4
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgcaccaccgaaacgactccagagctggggaagctgctgccgagagccgtgtgcaaatggcaataactattctctccaccttaccctgaacactgcagaccttggcggggaatgcatgaacaagctcaaatgcccgcctgctgaagcagagataatgaatcc
cccagcgcgtgggacttttccgctatcccggcittaggggatctcattaccctccaggggtcatctcatgacagccctcactccccgcagcagcctcagcagcctgacagacagtcggttcaaatgccaatcggcagactgccgagaatcattctctctctctctct
cctcaggggagctggcgagccaaccagccaagaagaggaaaaggtgtgggtctgctgcccctgcaagaggtcatcaactgtggcctgagcaggtgcaggaaccgaaaacgggacaccagatgcaaattagaaaatgtgaagagctaaagaaaaaacctggc
actcactagagagaacacctgttcccagcctgaagcattccgatggtcttttaa