

C18orf21-chromosome 18 open reading frame 21 Gene

Reference: PTXBC119766

| | | | |
|----------------|---|---------------------|----------------|
| Product type: | DNA | Storage Condition: | -20°C |
| Size/Quantity: | 2 ug | Delivery Condition: | +4°C, Blue Ice |
| Synonyms | UPF0711 protein C18orf21; HsT3108; PNAS-124; PNAS-131; XTP13; HBV X-transactivated gene 13 protein; HBV X-transactivated protein 13; HBV XAg-transactivated protein 13; chromosome 18 open reading frame 21 | | |

For research use only.

| | |
|-------------------------|---|
| GENE DESCRIPTION | chromosome 18 open reading frame 21 |
| GENE ID | 83608 |
| ORIGIN SPECIES | Human |
| ACCESSIONS | BC119766 |
| VECTOR | pDONR223 |
| NCBI SYMBOL | C18orf21 |
| LENGTH | ProteoGenix bp |
| SEQUENCE PRIMERS | Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG |

ORF (Only for reference, product provided without STOP codon)

atgagacagaagcactacctgaggtgcagcgcgggactgcacgacagctgccggccaagcccgtacctctctggcctacactctgcgcagatgataagagcacttttgaagaaacgtgtccactgtttccagctgttcttgataactctcagtgctcgtcacaacc
caaagccaggttgacaccacaaatacagaaactcttaatcgagaagcagagaactatacactcagtttaagaagcaaaaatggtgaaaaagttcaaaagctcaaaaagtgtattgtgatcactgttaaaatgcaacagaaacgtgaaacatcaggtaaaagtgaagctttgttca
acattgaagagcaatcctgccactcctacaagtaaacagcctgaagacaccagagagaaggaactgcaaacccaatcatgacatgtctggctcgaagagcaagagccagccatcggttttcagaacacactacatctggacagtcagatctactgtcctcaaaacaccagcaaaa
caaagaacacttctcaactaaaaatgttacttagtcagaatgaatccaaaagattcaaggtgacttcagaaatttctatctctctgaaggtggacttttaaaataa