

FLJ44881-FLJ44881 Gene

Reference: PTXBC119776

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms			
For research use only.			

GENE DESCRIPTION	FLJ44881
GENE ID	400661
ORIGIN SPECIES	Human
ACCESSIONS	BC119776
VECTOR	pDONR223
NCBI SYMBOL	FLJ44881
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgggtgctctgaattccatttccttagaaccacaagagctccttcagtgatgctggacaccttcagttgagcacagctcaggcaccgaaggacactccgctcatgtaagttggccttagtaggttttttccatgagtgaaatttcattcagagcaagtaaaacacagcctgtgacactgtt
tgtcccgaatgtaacagcgcctctctgtgatcagccgtgtacgagggtgtggaatacatgaacaacgaaacggaatatctgctctcatttctgtctttaatcattggctgctccactccagttacacctgggctgctctgatctgcagggttaactccaggggaacttcacctgtttgat
gacagccgtcaagtga