

C6orf122-chromosome 6 open reading frame 122 Gene

Reference: PTXBC119778

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	C6orf122; NCRNA00242; dJ266L20.5; long intergenic non-protein coding RNA 242		
For research use only.			

GENE DESCRIPTION	chromosome 6 open reading frame 122
GENE ID	401288
ORIGIN SPECIES	Human
ACCESSIONS	BC119778
VECTOR	pDONR223
NCBI SYMBOL	C6orf122
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgatccccggcactctcctcactcaggatcactcctcgccgcggttccagacgaggttccctcacggagcgagggaagattcaggcgtgtctgttctgtctttgtcagcacggacaaggatggagacccactggcagggcagccccgggtccctgttctcactt
caccactgggggaagcctgatccccatcgatcgcaaaggaacaaggaaagacacgtttacctggatggatgggtccctcacgggtggttggaaacccgttcagcgggagaggcagcgctcactcacgggtccccctggcgtgttaccaggagcgccacccgtgccc
cgccagagggaagtgcaggcgccagcacagcagggccgaaagccacactgtgggacccgcagcgagctccaagaaatccaagtcattcacaagaggagcttctccgcaaatcgttgaaaaataaacgagaaatgagtcgccaccggtatccgcttagtcagccgtacgaaa
acgcaacccggtcagctgacttctgctgccagccctcctcagtcagcaccggtcgcagaagggaaggcaggttag