

**INCA1-inhibitor of CDK interacting with cyclin A1 Gene**

Reference: PTXBC119781

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	protein INCA1; HSD45; inhibitor of CDK, cyclin A1 interacting protein 1		

**For research use only.**

---

<b>GENE DESCRIPTION</b>	inhibitor of CDK interacting with cyclin A1
<b>GENE ID</b>	388324
<b>ORIGIN SPECIES</b>	Human
<b>ACCESSIONS</b>	BC119781
<b>VECTOR</b>	pDONR223
<b>NCBI SYMBOL</b>	INCA1
<b>LENGTH</b>	ProteoGenix bp
<b>SEQUENCE PRIMERS</b>	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

**ORF (Only for reference, product provided without STOP codon )**

atgcaggtcaggatgatggagtaacatccccccttggcaagtgtccagggtgtcagccgatccaccccaagggtgcctccagagctcagaccatgccccagcgttatggagatgtctctggaagaacctaatcaaggccccccccactggctggagagcagc  
acattccaccatgctgagagccactgggtgtccagcttggctgtatctctctgagcagctccacccccgaaatgcttggagaagaagaaggagccatgtttgaaaggaatgacagcagcagggccttggggagctcccggcgggtgaggctgctactaccactgg  
aggacctaaagagcgtcagagcatcaacgaactgaagaagccagctggggcagctctgggctgcatctgagccagctgcttggcgaagaggctgtgatlccccagcaccatgaataccctgatctggaagagagagagcaacctatccacaggaaggaccgtt  
ttctactctggcaggccccagctgttggctcctggagccccctggatcaggaggagcttgcctcaggcagctgactctggcctgttcagcactgtcacagccagaaggaacccccctcacaatccctggggatggagttgcagcgtctgaagagtaa