

C2orf50-chromosome 2 open reading frame 50 Gene

Reference: PTXBC120930

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	uncharacterized protein C2orf50; chromosome 2 open reading frame 50		
For research use only.			

GENE DESCRIPTION	chromosome 2 open reading frame 50
GENE ID	130813
ORIGIN SPECIES	Human
ACCESSIONS	BC120930
VECTOR	pDONR223
NCBI SYMBOL	C2orf50
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atggggagccacccacccctgggtccagagaaccacatcagctgggtaccgattgccccccaccagccaccagcctcgggtctccccagctgccccgggtggccctatggccagcaggggtcttctggtggtggtccagggccccccaggtctgaaggcgagcgggtggtca
gggagctgcctgcgatggcgtgcagcagaccagctgtgggggagctcctggaggccgagcggcgggcccagcagcgtggattcagaactggagttcctgaaagactatgacccaatgggcaacaagaaggagcctgagaagttgccagaccatgtcctctcttgcggacac
agttccagttccacgaaccaggttggggcagcaggtggacacaccctgggacagactctcatccgcatggactctcttcacagaaggggccggagaagaagctggaggaccagatgcagccatctag