

LOC729233-FLJ40296 protein family member Gene

Reference: PTXBC120931

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	proline-rich protein 20B; FLJ40296 protein family member; proline rich 20B		
For research use only.			

GENE DESCRIPTION	FLJ40296 protein family member
GENE ID	729233
ORIGIN SPECIES	Human
ACCESSIONS	BC120931
VECTOR	pDONR223
NCBI SYMBOL	LOC729233
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atggaggaaccaaggccttcgaagcgaactcgtccatggccctaatcaagcctcaggtggcctctccagagccagcgtgtgttgcggaccctgaagcctcgtggaagcagatgggcccgcacagcccaacccgcaaaaccatcgcttacgtgaaaccttcagac
ggcagccccagctcgcagagagtcacccctctcgcagagagagccggcggggaggaagccggcggccagggcgagccgtggcagaaggcctgggccccgggggacgtggccagagacagggggcagaaggcttgatggcaccggacgtgcacatccaaact
ggaccacatggagagccagggccaccagggggaaccggaatcacggagaccgcagccttctcctttctgaaacaggtcctccgctggaactgtgcaggaaggccctggccccgacgtggcgcaacctgagctggggttcaggagccgcccgtgtcctgggctcaggctg
ttgactggcaaccgtcttgaccctctatccctgcacgggttttagggctcgtgggtgactcagcgggtttacaagtcattcaaacccccagggcacctacgtgcaagggtgccagtgttcctaccgacattgcgtattga