

TBX20-T-box 20 Gene

Reference: PTXBC120946

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	T-box transcription factor TBX20; ASD4; T-box protein 20; T-box 20		
For research use only.			

GENE DESCRIPTION	T-box 20
GENE ID	57057
ORIGIN SPECIES	Human
ACCESSIONS	BC120946
VECTOR	pDONR223
NCBI SYMBOL	TBX20
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atggagttcaggcgtccccaagcccaactctccgggccaacgcttcattgccgcgtcatgtcagcggcgctctaaaggagaaggagcgacggagaacacaatcaaacccctggagcaattgtggagaagtcgtcctgtgccagccctgggtgagctgacca
gcctggatgctcatgggagtttggaggcagtgccagcagcccgctcctcctctctgtgactgagccactgatccccaccaccccatcctccagtgaggaaatggccaaattgctgcagcctggagaccaaggagcttgggacaaattccatgagctgggcaccgaga
tgatcatcaccaagtcgggcaggaggtgtttccaaccatccgggtgctctttcgggggtggtcctgagggcgaagtacatagtcctgatggacatcgtccctgtggacaacaagaggtaccgctacgcctaccaccggctcctcctggctggtgctggcgaaggccgaccgcctgtgcc
agccaggtctctatgtgcatcagattctctttaccgggtgagcaactactcaaacagatgggtgtctttgaaaagggtgaaactcaccaacaatgaactggatcaaatggccataaatttgaactcaatgcataagtaccagccaagggtgcacatcattaagaagaagaccacacagcct
cattgctcaacctgaagctcgaagaatttagaactttcatcttccagaacaagtgtttacggcagtcactgcctaccagaatcaactgataacgaagctgaaaatagatagcaatcctttgccaaaggattccgggattcctcaggtcactgacattgagaggttaa