

GPR141-G protein-coupled receptor 141 Gene

Reference: PTXBC120951

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	PGR13; G-protein coupled receptor PGR13; G protein-coupled receptor 141		
For research use only.			

GENE DESCRIPTION	G protein-coupled receptor 141
GENE ID	353345
ORIGIN SPECIES	Human
ACCESSIONS	BC120951
VECTOR	pDONR223
NCBI SYMBOL	GPR141
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgcctggccacaatacctccaggaattcctcttgcgatccatagtgacaccccaatcaatcagcctctactcatagtgcttattggcgggctgggtgtcatttccattctttctcctgggtaaaatgaacacccggtcagtgaccaccatggcggtcattaacttgggtgtgtccacag
cggttttctgctgacagtgccatttgccttgacctacatcaagaagacttggatgtttgggtgcccctctgcaaaatttggagtgccatgctgcacatccacatgtacctcaggttctattctatgtgtgatcctgggtcaccagatccctcttcaagtgcagaaacaaatgggaattct
acagaaaactgcatgctgtggctgccagtgctggcatgtggacgctggtgattgtcatttgggtaccctggttctccgggtatggaatccatgaggaatacaatgaggagcactgttttaatttcacaaagagcttctacacatatgtgaaaatcatcaactatgtagtcatgtttgtca
tagccgttgcgtgattctgttggcttccaggttctcattatgttggatgagcagaagctacgccactcttactatcccaggagttctgggtcagctgaaaaacatttttatagggtcatccttgtttgttcttccctaccagttcttaggctattactgaatgtgtgacgcatc
caatgcctgtaacagcaagggttgcattttataacgaaatcttctgagtgtaacagcaattagctgctatgattgtcttcttcttcttgggggaagccattgtttaagcaaaagataattggcttatggaattgttttggccgttag