



Reference: PTXBC120959

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	KAP10.1; KAP18-1; KAP18.1; KRTAP10.1; KRTAP18-1; KRTAP18.1; keratin-associated protein 10-1; high sulfur keratin-associated protein 10.1; keratin-associated protein 10.1; keratin-associated protein 18-1; keratin-associated protein 18.1; keratin associated		

GENE DESCRIPTION	keratin associated protein 10-1
GENE ID	386677
ORIGIN SPECIES	Human
ACCESSIONS	BC120959
VECTOR	pDONR223
NCBI SYMBOL	KRTAP10-1
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

atgacgcgcctaccatctgccctgtctccagcgtctccgcactctggcaggctggaagctgcccagagagctgctgtagagcccccactgctgccctcagctgctgcccctfagccctgctgctgaccccagtagagccggtgctgctccagcccctgctgccaggcagcctgtgagccagccctgcacatcaggctgcaccagctctctgcacgccctgctgctgcagcagctctagctgcagcagccggctgctgcactctccccctgccagcaggcctgctgcatgcccgtctgctgcaagcctgtgtgctgctgctgccacctgctctaaaggattctctctatgctgcagcagctctagctgccagccaaactgctgtgcccctctctccagccagcagctctgctgtgctgctgttctgccaagcctgtgtgctatgtgcccaactgctctgagtcactctctctatgctgcagcagctctagctgccaatcagctgctgctgacctctccccctgccagcaggcctgctgctgcccctgctgcaagcctctgctgtgcccctctgctgtgcccctctgctgctggcctctccattctatgctgcagcagctctagctgctgacccactctgctgctgagacctctctctctgctgctctctctgcccctgtatgccagcccgctgctgtagtctctctctctgctgctccctgcctccagcagctgctgcccgcctctgctgctccctctctgcccgcctgctgctgctccctctctgcccgcctgctgctgcccctgctgctgcccgcctgctg