

C1orf77-chromosome 1 open reading frame 77 Gene

Reference: PTXBC120961

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	C1orf77; FL-SRAG; FOP; SRAG; SRAG-3; SRAG-5; pp7704; chromatin target of PRMT1 protein; friend of PRMT1 protein; small protein rich in arginine and glycine; chromatin target of PRMT1		

For research use only.

GENE DESCRIPTION	chromosome 1 open reading frame 77
GENE ID	26097
ORIGIN SPECIES	Human
ACCESSIONS	BC120961
VECTOR	pDONR223
NCBI SYMBOL	C1orf77
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atggctgcacagtcagcgcgaaagtgtgtctaaaaagcaccaccaagatgtctctaatgagcgttactaatatgtctgaagaacaacagccgacgccagtgaatttcgggcttcgatgcagcaacaacagcagctagccagtgccagaacagaagactggccagcagatgga
gaatagacctctgtccaggcagcattaaaacttaagcagaagagcttaagcagcgcctgggtaagagtaacatccaggcacgggttaggcgacccatagggccctggccaggggagcaatcggaggacgaggcctaccataatccagagagcctggccagaggaggactac
gtgggggacgtgccaccagaaccctacttagggcggtatgtcactccgaggtcaaaacctgtctccgaggtggacgagccgtagctccccgaatgggcttaagaagaggtgtgttcgaggtcgtggaggtcctgggagagggggcctagggcgtggagctatgggtgtggcgga
atcgggtgtagaggtcgggtatgataggtcggggaagaggggctttggagccgagggccgagccgtggacgagggagaggtgccctgtctgcctgtattgaccaaggagcagctggacaaccaattggatgcataatgtcgaacaaaaggacacctggatgtgagttg
gatgcctacatggcgagacagatcccgaaccaatgattga