

FBXO32-F-box protein 32 Gene

Reference: PTXBC120963

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	Fbx32; MAFbx; F-box only protein 32; atrogin 1; muscle atrophy F-box protein; F-box protein 32		

For research use only.

GENE DESCRIPTION	F-box protein 32
GENE ID	114907
ORIGIN SPECIES	Human
ACCESSIONS	BC120963
VECTOR	pDONR223
NCBI SYMBOL	FBXO32
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgaatatttggaaaaagtgtactgaaagtccttgaagaccagcaaaacattagactaataagggaactactcagaccctctacacatccttatgtacactggtccaaagagtcggcaagtcgtgctggtcgggaacattaacatgtgggtgtatcggatggagacgatttccactggcagcagcgtgaacaacattcagatcaccaggcctgcttcaaaaggcctcactcactgacctgcttctgctcacaactgaacatcatgcagaggctgagcgacggcgggacgtggtcagcctgggcccagctgccccgacctgcacgtgctcagcgaagaccggctgctgtggaagaaactctgccagtaccacttctccgagcggcagatccgaaacgattaattctgtcagacaaaggcagctggattggaagaagatgtatttcaaaactgtccgatgttacccaaggaaagagcagatggagatacccttcagctctgcaaacactgtcacatccttctggaaggga