

LOC57228-small trans-membrane and glycosylated protein Gene

Reference: PTXBC120965

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	hSMAGP; small cell adhesion glycoprotein; small cell transmembrane and glycosylated protein; small trans-membrane and glycosylated protein; small transmembrane and glycosylated protein		

For research use only.

GENE DESCRIPTION	small trans-membrane and glycosylated protein
GENE ID	57228
ORIGIN SPECIES	Human
ACCESSIONS	BC120965
VECTOR	pDONR223
NCBI SYMBOL	LOC57228
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgaccagcctcctgactactccttctcaagagaagaactgatgaccacccaattttacagccactgaggccctgtcccccagaagatggagccagcagcactcattgcagttgttatcaccgtgtcttctcaccctgctctcggtctgactctgatctctttacgtgtacaagaaca
aaggcagctactcacctatgaacctacagaaggtgagccagtgccatcgtccagatggagagtgacttggccaagggcagcgagaaagaggaaatattcatctaa