



Reference: PTXBC120975

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	MUC1/ZD; ADMCKD; ADMCKD1; CA 15-3; CD227; EMA; H23AG; KL-6; MAM6; MCD; MCKD; MCKD1; MUC-1; MUC-1/SEC; MUC-1/X; PEM; PEMT; PUM; mucin-1; H23 antigen; breast carcinoma-associated antigen DF3; cancer antigen 15-3; carcinoma-associated mucin; episialin; krebs		

For research use only.

GENE DESCRIPTION	mucin 1, cell surface associated
GENE ID	4582
ORIGIN SPECIES	Human
ACCESSIONS	BC120975
VECTOR	pDONR223
NCBI SYMBOL	MUC1
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

atgacacggcgaccacgctctcttctctctgctgctgctctccacagtcttacacgtaccacagcccctaaccgcacaagtgttacgggtctggctatccaagctctaccacgggtggagaaaaggagacttcgctaccagagaagttcagtgcccagctctactgagaagaat
gctttaattctctctctggaagatcccagaccgactactaccaagagctcgagagagacattctgaaatgttttcgagattataaacaagggggtttctggccctcctaataatagttcagccagagatctgtgtggtaacaattgacttcggcctccgagaaggtaccatcaatgtcca
cgactgtggagacacagttcaatcagfataaaacgggaagcagctctcgtatataacctgacgatctcagacgtcagcgtgagtgatgtgccattctctctcgtccagctctggggctggggctggcagcctggggcattcgcgctgctgctgctgtgttctgggtgcgctggggcattgtc
tatctcattgccttgctgctgctgacgtcccgccgaagaactacggcagctcgacatctttccagccgggatacctaccatctcagcgcagtagtaccacacccatccacacccatggcgctatgtgccccctagcagtaccgatcgtagccctatgagaaggtttctcaggttaattgtt