

NXT2-nuclear transport factor 2-like export factor 2 Gene

Reference: PTXBC120984

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	P15-2; NTF2-related export protein 2; protein p15-2; nuclear transport factor 2 like export factor 2		

For research use only.

GENE DESCRIPTION	nuclear transport factor 2-like export factor 2
GENE ID	55916
ORIGIN SPECIES	Human
ACCESSIONS	BC120984
VECTOR	pDONR223
NCBI SYMBOL	NXT2
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgagaaaatacagaagccactggtctcaggagacagagacgataccaaagaagaagtaactattatgagggccacacaccagccactcaagccctgcggaccggacacgtgaggaggtagtgaacggacactgccagaacacactgctacaagggtccagatggccacgtc
tctggattttaaacttatgtatcagggcatgtagagctgctgaggagttgtcaatattactatgagacaatggataaaagaagacgggcactaaccaggtgtatctggacaagggccacttaatatggaatggaatgctgtttcagggtggtatgccctaaataatttttgacacattgc
cttctagtgtggtccaggtcaatatgttagattgccaaccaggtcatgagcaagcaactcagtcaccaactacagttctgtgtgaccagtggaactgtgaagtttgatggaaacaacaacatttctcaaccagaactctgctgactgtcagtcactcccaacaatactgtgtggaagatt
gcaagtgattgcttccgttticaagattggtctagtagtaa