

RAB6C-RAB6C, member RAS oncogene family Gene

Reference: PTXBC120999

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	RAB6C, member RAS oncogene family; ras-related protein Rab-6C; rab6-like protein WTH3		

For research use only.

GENE DESCRIPTION	RAB6C, member RAS oncogene family
GENE ID	84084
ORIGIN SPECIES	Human
ACCESSIONS	BC120999
VECTOR	pDONR223
NCBI SYMBOL	RAB6C
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgtccggggcgggagacttcgggaatccgctgaggaaattcaagctgggtgttcctgggggagcaagcgttgcagaagacatctttgatcaccagattcaggtatgacagttttgacaacacctatcaggcaataattgcattgacttttatcaaaactatgtacttggaggatggaacaat
cgggcttcggctgtgggatacggcgggtcagggaacgtctccgtagcctcattccagggtacatccgtgattctgctgcagctgtagtttacgafatcacaaatgtaactcattcagcaaaactacaaagtggattgatgatgcagaacagaaagggaagtgtatcatcacgctag
taggaaatagaacagatcttctgacaagaggcaagtgtcagttgaggaggagagaggaaagcgaaggctgaatgttacgtttattgaaactaggcgaataatgtaaacgagctcttgcagctgtgtagcagcagctttgccgggaatggaagcacacagcagcga
agcagagaagacatgagtgacataaaactggaaaagcctcaggagcaaacagtcagcgaagggtgttcctgctactctccatgcatctcaaccctcctcagaagccccctactcttcattgactgcagtgatgaattggcctgaaccttttccttcattaataacgttttcaattca
tcattgctgcctgtctcgtggaggtga