

C21orf58-chromosome 21 open reading frame 58 Gene

Reference: PTXBC121008

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	uncharacterized protein C21orf58; chromosome 21 open reading frame 58		

For research use only.

GENE DESCRIPTION	chromosome 21 open reading frame 58
GENE ID	54058
ORIGIN SPECIES	Human
ACCESSIONS	BC121008
VECTOR	pDONR223
NCBI SYMBOL	C21orf58
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atggcgcgatctcggctccctgcaacctccctcaggaagccgtggaagctcgaccgcagaaactctctctgaactcagccacagctctctgtgtgctggtctccaggagtaaggcccgccctgcaggcaaacccggtgcttggctcctctgagcagttcttctcgagtaa
cagaaccaggaggagggtggctgtgctccctgccccctacagctctctctcagctcccacatgctggactcagcagcagcaagtgaaccgactgacgtgaagctctggacagaagctggagcaagaacggcagaacgtggaaggggacctgagggcctc
cacctcagccaggaatgaggaccggccggacgatccctgcagactgcttgaagagaaggaggacctctgcagagactcgggtaggccctgcacccttgagcgaacctggctgagctcactgtag