



Reference: PTXBC121015

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	ZNT1; ZRC1; zinc transporter 1; solute carrier family 30 (zinc transporter), member 1; znT-1; solute carrier family 30 member 1		

For research use only.

GENE DESCRIPTION	solute carrier family 30 (zinc transporter), member 1
GENE ID	7779
ORIGIN SPECIES	Human
ACCESSIONS	BC121015
VECTOR	pDONR223
NCBI SYMBOL	SLC30A1
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

[illegible]

More info: Please contact us for any question: contact@proteogenix.fr