



Reference: PTXBC121017

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	UPF3B pseudogene 3; UPF3B pseudogene 2; UPF3B pseudogene 1; HUPF3B; MRX62; MRXS14; RENT3B; UPF3BP1; UPF3BP2; UPF3BP3; UPF3X; Upf3p-X; regulator of nonsense transcripts 3B; mental retardation, X-linked 62; nonsense mRNA reducing factor 3B; up-frameshift su		

For research use only.

GENE DESCRIPTION	UPF3 regulator of nonsense transcripts homolog B (yeast)
GENE ID	65109
ORIGIN SPECIES	Human
ACCESSIONS	BC121017
VECTOR	pDONR223
NCBI SYMBOL	UPF3B
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

[illegible]