

WFDC10B-WAP four-disulfide core domain 10B Gene

Reference: PTXBC121071

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	protein WFDC10B; WAP12; protease inhibitor WAP12; WAP four-disulfide core domain 10B		

For research use only.

GENE DESCRIPTION	WAP four-disulfide core domain 10B
GENE ID	280664
ORIGIN SPECIES	Human
ACCESSIONS	BC121071
VECTOR	pDONR223
NCBI SYMBOL	WFDC10B
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atggcaccacagactctgctgctgtctgtctgtctgtctgtctgtcgcaggccaggaggataccgtgacaagatgaggatgcagagaatcaaggctctgtgagaagcgaccagcatagatctatgcattccaccactgttcatttccaaaagtgtgaacaataagatatgctgttcagcctctgtgggaacatttgtatgagcatcctatga