

LOC729862-similar to Striatin Gene

Reference: PTXBC121074

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms			
For research use only.			

GENE DESCRIPTION	similar to Striatin
GENE ID	729862
ORIGIN SPECIES	Human
ACCESSIONS	BC121074
VECTOR	pDONR223
NCBI SYMBOL	LOC729862
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atggccgtggccagtacgaagagtcggggagacgggtgaggtacaggctcagctgcggccaagactctgctcgaaggatattgtggctggagacatgagcaacaaagcttctgggagcagaaggagggtccaagacctcatcaacaattaaggccagattgccttcctaca
gggagaagggaaggccaagaaaattgaagaaggatctgtgaggatgatcaaatgttgagatgctcttaaacagaaaagaccaaataccacaagttgaaatatgggacagaattgaatcaggagctatgaagcctccaagctatgattctgatgaagtaataagaagtg
cagccacaacaaacagccagttaatgtga