

APOC3-apolipoprotein C-III Gene

Reference: PTXBC121081

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	APOCIII; HALP2; apolipoprotein C-III; apo-CIII; apoC-III; apolipoprotein C3		

For research use only.

GENE DESCRIPTION	apolipoprotein C-III
GENE ID	345
ORIGIN SPECIES	Human
ACCESSIONS	BC121081
VECTOR	pDONR223
NCBI SYMBOL	APOC3
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgcagccccgggtactcctgtgtgtgcctcctggcgctcctggcctctgcccagcttcagagcgaggatgcctccctctcagcttcacatgagggtacatgaagcacgccaccaagaccgccaaggatgactgagcagcgtgcaggagtcaccaggtggcccagcaggccag
gggctgggtgaccgatggcttcagttccctgaagactactggagcaccgttaaggacaagtctctgagttctgggattggaccctgaggtcagaccaacttcagccgtggctgcctga