

LOC151300-hypothetical LOC151300 Gene

Reference: PTXBC121084

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms			
For research use only.			

GENE DESCRIPTION	hypothetical LOC151300
GENE ID	151300
ORIGIN SPECIES	Human
ACCESSIONS	BC121084
VECTOR	pDONR223
NCBI SYMBOL	LOC151300
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgttactggaggagataggctgtgtcatcctgcatgggaatgacctaatacaaaagtatcagtgactgctcaggactgtgaaaggtcacttggagcagtgaggaaaggtctggcagaggaataagaagatacaggcttagaacaagccctgggactcaaccaagtcattaac
atctggggccttgcgtttattatcctccagcagataaagcaataacctgccctgtcttctaaccttggaggggtcaaaaaagtatgggtgtgaaagcattttgtaactgtag