

LOC283332-hypothetical protein LOC283332 Gene

Reference: PTXBC121096

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	uncharacterized LOC283332		
For research use only.			

GENE DESCRIPTION	hypothetical protein LOC283332
GENE ID	283332
ORIGIN SPECIES	Human
ACCESSIONS	BC121096
VECTOR	pDONR223
NCBI SYMBOL	LOC283332
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgcatggcaacagaccaggcctgtgccaggcctaccaaggagtgagctgacggctgagagttttggaatgtgaggtgggatgacaagcctggggcctgtcttctgtctcagtgaggactgtccagcctcagcccaactccagcaaaagccaggtgtggggagca
ggggagggtgagacagagatgaggaggacagaaatgtgtttcacctgacccagcaattcctgcagctgggaagcatggcgactcagaaagccgatgcctgcccagagcttcggcttctgtcccgcgttttggagaccatctctgagagggcacagtgta