

FGF6-fibroblast growth factor 6 Gene

Reference: PTXBC121097

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	HBGF-6; HST2; fibroblast growth factor 6; FGF-6; HST-2; HSTF-2; heparin secretory-transforming protein 2; heparin-binding growth factor 6		

For research use only.

GENE DESCRIPTION	fibroblast growth factor 6
GENE ID	2251
ORIGIN SPECIES	Human
ACCESSIONS	BC121097
VECTOR	pDONR223
NCBI SYMBOL	FGF6
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atggccctgggacagaaactgttcatcatatgtccggggagcaggacgtctgcagggcacgtgtgggtctctgtcttcttaggcacatctagtgggcatgtgtgtccctgcctgcaggcaccgtgccacaacacgtctggactcaggggtggggcacctgtgtccag
gtctcgcgcgggtagctggagagattgccgggtgaactgggaaagtggctatttgggtgggatcaagcggcagcggggtctctactgcaacgtgggcatcggctttcacctccaggtgtctccccgacggccggtatcagcgggaacccacgaggagaacccctacagctgtctgg
aaatttccaccgtggagcagggcgtggtgagtcctttggagtgagaagtgcctctctgtcgccatgaacagtaaaaggaagattgtacgcaacgccagcttccaagaagaatgcaagttcagagaacccctctgcccacaattacaatgcctacgagtcagactgtaccaaggac
ctacattgccctgagcaaatacggacgggtaagcggggcagcaaggtgtccccgatcatgactgtcactcatttcttcccaggatctaa