

MGC3771-hypothetical LOC81854 Gene

Reference: PTXBC121099

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms			
For research use only.			

GENE DESCRIPTION	hypothetical LOC81854
GENE ID	81854
ORIGIN SPECIES	Human
ACCESSIONS	BC121099
VECTOR	pDONR223
NCBI SYMBOL	MGC3771
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atggggtttctccatgtcggtcaggctggctcgaactcctgacctcaggtgatccgccacctcagcctcctaa