

HOXD12-homeobox D12 Gene

Reference: PTXBC121103

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	HOX4H; homeobox protein Hox-D12; Hox-4.7, mouse, homolog of; homeo box D12; homeobox protein Hox-4H; homeobox D12		

For research use only.

GENE DESCRIPTION	homeobox D12
GENE ID	3238
ORIGIN SPECIES	Human
ACCESSIONS	BC121103
VECTOR	pDONR223
NCBI SYMBOL	HOXD12
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgtgtgagcgcagctctacagagcgggtatgtggctcgtcttgaatctgcagtcgcagactcttctacttctccaactgagccgaatggcgccagttggcgcgcttccccctatctctaccgcggcgcgcgtgccctgggcccacagcccgctctgcgccccg
cgagcctgcggcgccactgccttgcggcgtcttcgcagccctacctggctcggcctctcggcctgcagcccccaacagccaaagacggaccggaagagcaggctaaagtctatgcgccgaagcggcgtggccagaggagcgggtgtacctggcgtcttc
gccccgagctctagcctggctcgtcagtggtcgtctcaaaagcggccaagtatgactacgctggtgtgggtcgtgccacgccgggtccacgacctgctccaggggctccctgcgccctggcttcaaggacgacaccaaggcccgctcaactgaacatgacagtgcaggcgg
cggcggttgcctcttgcctgcgaccttactgccgacggtgtaaacggtgcccatgctccccggccggtttggcgccggatgggaggtgggttcaaggagagtgtaaggagggtga